

TITRATOR

Mettler Toledo T70 Titrator, Velika Britanija



Lokacija: Katedra za farmacevtsko tehnologijo

Titratore omogoča varno, prilagodljivo in učinkovito izvajanje titracij ob simultanjem mešanju. Aparatura nam nudi tudi izvajanje titracij ob vzdrževanju pH (pH-Stat titracije) ter s stalnim nadzorom temperature.

Pogosta uporaba v farmaciji:

- *In vitro* lipoliza
- Vsebnost prostih maščobnih kislin v oljih
- Določanje trdote vode
- Merjenje kislinkega števila
- Vsebnost vitamina C
- Titracija površinsko aktivnih snovi
- Idr.

TITRATOR

Mettler Toledo T70 Titrator, UK



Location: Department of Pharmaceutical Technology

The titrator enables safe, adaptable and efficient performance of titration with simultaneous stirring. The machine offers also the application of pH-Statting titration as well as continuous temperature control.

Scope of the titrator usage in Pharmacy:

- *In vitro* lipolysis
- Free fatty acid content in oils
- Acid number determination
- Total hardness of water
- Vitamin C content
- Surfactant titration
- Others