

GLATT GMPC I

Glatt GmbH, Binzen, Nemčija



Naprava se uporablja za oblaganje tablet.

Uporaba:

- možnost uporabe organskih topil
- boben velikosti 0,8 ter 1,6 L
- možnost zajema podatkov tekom oblaganja
- možnost oblaganja tablet s premeri nad 3 mm

Lokacija: Katedra za farmacevtsko tehnologijo

GLATT GMPC I

Glatt GmbH, Binzen, Nemčija



The device is used for tablet coating.

Applications:

- Vessel size: 0.8 and 1.6 L
- Process data logging
- Possible to coat tablets with diameter > 3 mm
- Possible to coat using organic solutions

Location: Department of Pharmaceutical Technology