

Tabletirka na udarec

IMA Killian SP300



Lokacija:

Katedra za farmacevtsko tehnologijo

Tabletirka na udarec SP300 je popolnoma instrumentirana in omogoča spremljanje sile na zgornjem in spodnjem pečatu, sile izmeta in odmika zgornjega pečata. Primerna je za napreden študij kompresibilnosti in analize krivulj sila od odmika. Opremljena je s polnilnim čolničkom na gravitacijsko in prisilno polnjenje. Omogoča širok spekter aplikacije od ročne izdelave le nekaj tablet do izdelave serij polindustrijskih velikosti.

Odljuje jo:

- popolna instrumentiranost
- visoka fleksibilnost

Na tabletirko je mogoče namestiti najpogosteje uporabljene EU-B in EU-D pečate.

Single-punch tableting press

IMA Killian SP300



Single-punch tableting press SP300 is fully instrumented machine and enables measurements of upper and lower punch forces, ejection force and displacement of the upper punch. It is well suited for advanced studies of compressibility, as well as analysis of force-displacement curves. It is equipped with gravity and forced filling shoes. The machine has a wide spectrum of applicability from manual production of only a few tablets to batch production on semi-industrial scale.

Qualities:

- full instrumentation
- high flexibility

Most often used EU-B and EU-D tooling can be mounted on the tableting press.

Location:

Department of Pharmaceutical Technology